

# Spatially and temporally resolved fs/ps CARS measurements of rotation-vibration non-equilibrium in a CH<sub>4</sub>/N<sub>2</sub> nanosecond-pulsed discharge

TIMOTHY Y. CHEN<sup>1</sup>, BENJAMIN M. GOLDBERG<sup>2\*</sup>,  
CHRISTOPHER J. KLIEWER<sup>2</sup>, EGEMEN KOLEMEN<sup>1,3</sup>, AND YIGUANG JU<sup>1</sup>

---

AIAA SCITECH FORUM, JANUARY 11-21, 2021 (PDL-10)

<sup>1</sup> *Department of Mechanical and Aerospace Engineering, Princeton University, Princeton, NJ 08544*

<sup>2</sup> *Sandia National Laboratories, Livermore, CA 94550*


<sup>3</sup> *Princeton Plasma Physics Laboratory, Princeton, New Jersey, 08540, USA*

*\*Presently at Lawrence Livermore National Laboratory*




# Question: How does molecular energy transfer and methane conversion evolve in space and time in plasma methane reforming?

- Want to utilize plasma to electrify  $\text{CH}_4$  conversion technologies
- Need to understand energy transfer to and from molecules
  - How does rotation-vibrational non-equilibrium evolve in time and space? How does this influence  $\text{CH}_4$  conversion?
- Need *in situ* measurements of rotation-vibration non-equilibrium and molecular number density  $\rightarrow$  hybrid fs/ps CARS




# What is hybrid fs/ps CARS?

- Use broadband femtosecond pulse for simultaneous acquisition of many Raman transitions and picosecond probe for spectral resolution
- Several advantages:
  - $I_{CARS} \propto I_{fs}^2 * I_{ps}$  → ultrafast pulses provide high signal and straightforward spatially-resolved one-dimensional (1-D) imaging
  - No need to scan Stokes frequency and only needs 2 beams
  - Delayed probe pulse avoids non-resonant background
  - Can make near surface measurements (within 50 to 100  $\mu m$ )
- Would like to make spatially-resolved 1-D measurements in pulsed plasma



# First level vibrational temperature can be measured from the pure rotational CARS spectrum



Increasing vibrational energy level red-shifts the rotational energy level  $F(v,J)$


$$F(v,J) = B_v(J(J+1)) - D_v(J^2(J+1)^2)$$

$$B_v = B_e - \alpha\left(v + \frac{1}{2}\right) + \gamma_e\left(v + \frac{1}{2}\right)^2, D_v = D_e + \beta_e\left(v + \frac{1}{2}\right)$$

More details available in  
 Chen T.Y. et al. Opt. Lett. (2020)



# Rotational CARS can easily be extended to 1-D with only two beams for time and spatially resolved measurements





- No need for frequency-converting optics (i.e. OPA) for separate Stokes beam [1]
  - Bandwidth of 100 fs Ti:Sapphire regenerative amplifier enough to resolve rotations ( $<200 \text{ cm}^{-1}$ )




# Discharge geometry and experimental setup

- 40% CH<sub>4</sub>, 60% N<sub>2</sub>, 60 Torr
- 4 kV, 500 ns pulse width, 20 Hz, 8 mm gap, 220 Ω resistor
- fs laser: <7 fs, 0.6 mJ (bandwidth up to CH<sub>4</sub> ν<sub>1</sub> Q-branch at 2916 cm<sup>-1</sup>)
- ps laser: 65 ps, 7 mJ
- Want to measure rotational N<sub>2</sub> CARS and vibrational CH<sub>4</sub> CARS
- All data is time-resolved and 1-D with 40 μm resolution starting within 150 μm of cathode

Generates stable thin filament (d ~ 1mm)




# Significant spatial gradients in vibrational temperature of $N_2$ were detected as a function of time



1500 K set as detection limit for evaluation of  $T_{\text{vib}}$



V-T relaxation and peaks in  $T_{\text{vib}}$  are immediately obvious from time evolutions of  $T_{\text{vib}}$  and  $T_{\text{rot}}$




$T_{\text{vib}}$  is strongly peaked  $\sim 1$  mm from the cathode






There is some spatial structure to the CH<sub>4</sub> profile but smoother than N<sub>2</sub> T<sub>vib</sub> and T<sub>rot</sub>



$I_{\text{CARS}} \propto N^2$  so  $\sqrt{\frac{I}{I_0}}$  is a  
 measure of N or number  
 density




There is some spatial structure to the CH<sub>4</sub> profile but smoother than N<sub>2</sub> T<sub>vib</sub> and T<sub>rot</sub>




- Around 1 mm there is a local dip in CH<sub>4</sub> number density where peak in T<sub>vib</sub> was observed
- How does this correlate with temperature? ( $T \propto \sim 1/N$  for constant pressure ideal gas)




# Time evolution shows significant number density drop unrelated to temperature (high CH<sub>4</sub> conversion)



# CH<sub>4</sub> number density appears to be insensitive to high N<sub>2</sub> T<sub>vib</sub>



Time interval with high T<sub>vib</sub>





CH<sub>4</sub> number density is almost flat with near constant T<sub>rot</sub>



# Some remaining questions

- Why is there this highly peaked structure in  $T_{\text{vib}}$ ?
  - Possible indication of transition or formation of a pulsed spark and local non-uniformities in electron density
- Why is the peaked structure not more pronounced in the  $\text{CH}_4$  number density figure?



# Conclusions and future work

---

- Measured  $T_{\text{vib}}/T_{\text{rot}}$  and  $\text{CH}_4$  number density in a 40%/60%  $\text{CH}_4/\text{N}_2$  pin-to-pin pulsed discharge with time and spatially-resolved broadband fs/ps CARS
- Spatial structure of  $T_{\text{vib}}$  has a large peak and is non-uniform
- Majority of  $\text{CH}_4$  consumption likely due to electron-impact and electronically excited  $\text{N}_2$  quenching
- Peak in  $T_{\text{vib}}$  found to correlate with a dip in  $\text{CH}_4$  number density
- 2-D modeling of the discharge formation is underway to try to understand the spatial structure in  $T_{\text{vib}}$  and  $\text{CH}_4$  number density



# Acknowledgements

---

This material is based upon work supported by the U.S. Department of Energy (DOE), Office of Science, Office of Workforce Development for Teachers and Scientists, Office of Science Graduate Student Research (SCGSR) program. This work was also supported by ExxonMobil through its participation in the Princeton E-filiates Partnership of the Andlinger Center for Energy and the Environment. TC was partially supported through the Program in Plasma Science and Technology at Princeton University (PPST) Fellowship. CK and BG were supported by the Office of Chemical Sciences, Geosciences, and Biosciences, Office of Basic Energy Sciences, U.S. Department of Energy. Sandia National Laboratories is a multi-mission laboratory managed and operated by National Technology and Engineering Solutions of Sandia, LLC, a wholly owned subsidiary of Honeywell International, Inc., for the U.S. Department of Energy's National Nuclear Security Administration. Any subjective views or opinions that might be expressed in this presentation do not necessarily represent the views of the US DOE or the US Government.

We would like to thank Dr. Benjamin Yee for lending us the power supply used in this work.

## Questions?

